

DIRECTORY OF SERVICES

Area of analysis: *in vitro* evaluation of medical devices

Mode of Analysis:

Quality tests

Procedure	Sample material	Method	SOP/Version	Instrument
Batch Release of Heparin Adsorbers	Adsorber devices	Rezirculation model using heparinized human blood	AA-063/Version 01/09.2008	
Evaluation of Redness of Hemodialyzers after use	Hemodialyzers	Image analysis	AA-067/Version 01/02.2009	Minolta Chroma-Meter CR-310